



Concentrated Lemon Zip Disinfectant

Disinfectant, Cleaner, Mildewstat, Fungicide, Sanitizer, Virucide*, Deodorizer.

Description

This product is a no-rinse disinfectant cleaner that disinfects, cleans and deodorizes in one labor-saving step. This product is for use on hard, non-porous surfaces in institutional facilities, businesses, office buildings, public restrooms, hotels, motels, hospitals, and healthcare facilities. This product can be used on washable hard, non-porous surfaces such as: glass surfaces, metal, plastic, sealed stone, and vinyl.

Directions

FOR USE AS A ONE-STEP, GENERAL, HOSPITAL, MEDICAL DISINFECTANT, FUNGICIDE, VIRUCIDE*, CLEANER:

1. Pre-clean heavily soiled areas.
2. Apply use solution of 2 oz. of this product per gal. of water to disinfect hard, non-food contact, non-porous surfaces with a sponge, brush, cloth, mop, by immersion, mechanical spray device. For spray applications, spray 6-8 inches from surface. Do not breathe spray.
3. Treated surfaces must remain wet for 10 minutes.
4. Allow to air dry. Rinsing of floors is not necessary unless they are to be waxed or polished.
5. Prepare a fresh solution daily or when visibly dirty.

For detailed instructions and safety information, please review container label and SDS.



Features

- **Mold and mildew control**
- **Deodorizer**
- **Glass Cleaner**

Applications

This product can be used on washable hard, non-porous surfaces such as:

- Glass surfaces
- Metal, sealed stone
- Plastic, vinyl

Specifications

Color	Lemon
Fragrance.....	Yellow
pH Range	11.5 - 12.0
Dilution	2 oz./gal.
Density.....	8.34
EPA Registered	Yes

Packaging

ITEM #	SIZE
BSL65180065T	4/64 ounce



Contact your Account Representative for more information

140 Private Brand Way, Athens, TN 37303 | mybrightsolutions.com

In case of emergency call Infotrac 800.535.5053

© 2018 Midlab, Inc. All right reserved.

BSL6518 | 03192019 - ESL 070218